


MOONEY VISCOMETER REPORT

ABC

TEMP(C)	100.0
PREHEAT TIME(MIN)	1.0
TORQUE UNIT	Mooney



ID1	ID2	ID3	QC	min	max	ML(1+4)	Date
345			OK	82.7	95.1	82.7	2016/04/06 15:33:09
567			OK	80.1	111.1	90.3	2016/04/07 15:00:51
789			OK	82.4	95.8	87.7	2016/04/12 10:56:37
999			OK	88.6	109.3	88.6	2016/04/12 12:02:19

Lower Limit			80.50
Upper Limit			92.50
Average	83.4	102.8	87.32
Standard deviation	3.62	8.55	3.27
Coefficient of variation(%)	4.34	8.32	3.74