


# RHEOMETER REPORT ( SAMPLE : H )

TEMP(C)	177.0
TIME(MIN)	6.0
TORQUE UNIT	kgf*cm



試験機状態	測定No.	ML	MH	ts1	ts2	tc10	tc50	tc90
リビルト前	1	2.02	19.50	1.50	1.72	1.67	2.13	2.47
	2	2.05	19.74	1.48	1.70	1.67	2.13	2.45
	3	2.03	19.47	1.52	1.72	1.68	2.15	2.48
リビルト後	1	2.03	19.53	1.49	1.71	1.67	2.13	2.46
	2	2.01	19.29	1.53	1.74	1.69	2.16	2.49
	3	2.03	19.53	1.52	1.75	1.70	2.18	2.52