


JIS K 6251 Tensile



ID1/ID2/m ² (Area)	MPa(Max)	%(Max)	MPa(50%)	MPa(100%)	MPa(150%)	MPa(200%)	MPa(250%)	MPa(300%)	Date
123/C/10.00	5.70	440.95	1.19	2.11	2.89	3.53	4.04	4.49	2020/06/12 16:51:59
123/B/10.00	5.57	426.92	1.27	2.17	2.94	3.51	3.99	4.43	2020/06/12 16:53:13
123/A/10.00	5.70	444.70	1.19	2.03	2.78	3.40	3.92	4.38	2020/06/12 16:54:20

Lower Limit									
Upper Limit									

Median	5.70	440.95	1.19	2.11	2.89	3.51	3.99	4.43
Average	5.66	437.52	1.22	2.10	2.87	3.48	3.98	4.43
Standard Deviation	0.08	9.38	0.05	0.07	0.08	0.07	0.06	0.06
CV(%)	1.33	2.14	3.87	3.50	2.92	2.05	1.44	1.29